Instructions for Use - HAART 200 Aortic Annuloplasty Device

Instructions for Use - HAART 200 Aortic Annuloplasty Device	2
Инструкции за употреба – изделие за аортна анулопластика HAART 200	8
Upute za upotrebu – Proizvod HAART 200 za aortnu anuloplastiku	16
Návod k použití - Zařízení k aortální anuloplastice HAART 200	23
Brugsanvisning - HAART 200 aortisk annuloplastikanordning	29
Instructies voor gebruik - HAART 200-aorta-annuloplastieimplantaat	35
Käyttöohje - HAART 200 -aortta-annuloplastiarengas	42
Mode d'emploi - Dispositif d'annuloplastie aortique HAART 200	48
Gebrauchsanweisung - Aortenannuloplastie-Implantat HAART 200	55
Οδηγίες χρήσης - Συσκευή αορτικής δακτυλιοπλαστικής HAART 200	62
Használati utasítás - HAART 200 aortabillentyű anuloplasztikai eszköz	70
Istruzioni per l'uso del dispositivo per anuloplastica aortica HAART 200	77
Bruksanvisning - HAART 200 aortisk annuloplastikkenhet	84
Instrukcja użytkowania - Urządzenie do annuloplastyki aortalnej HAART 200	90
Instruções de utilização – Dispositivo de anuloplastia aórtica HAART 200	97
Instrucțiuni de utilizare – Dispozitiv de anuloplastie aortică HAART 200	104
Návod na použitie – Pomôcka na aortálnu anuloplastiku HAART 200	111
Navodila za uporabo – naprava za anuloplastiko aortne zaklopke HAART 200	117
Instrucciones de uso – Dispositivo de anuloplastia aórtica HAART 200	123
Bruksanvisning - HAART 200 aortaannuloplastikenhet	130
Kullanım Talimatları - HAART 200 Aort Anüloplastisi Cihazı	136


Instructions for Use - HAART 200 Aortic Annuloplasty Device

Table of Contents

1. Device Labeling Symbols	2
2. Intended Use	2
3. Indications for Use	2
4. Annuloplasty Device Description	2
4.1. Overview	2
4.2. Technological Characteristics	2
4.3. Accessories	3
5. Contraindications	3
6. Warnings	3
7. Precautions	3
8. Magnetic Resonance (MR) Safety	4
9. Potential Adverse Events	4
10. How Supplied	4
10.1. Packaging	4
10.2. Storage	4
11. Directions for Use	5
11.1. Sizing	5
11.2. Handling and Preparation Instructions	5
11.3. Device Implantation	5
11.4. Sterilization	7
Disclaimer of Warranties	7
Patents	7

1. Device Labeling Symbols

2. Intended Use

The HAART 200 Aortic Annuloplasty Device is intended to be used to correct annular dilatation and/or maintain annular geometry of the aortic valve in adult patients undergoing bicuspid aortic valve repair.

3. Indications for Use

The HAART 200 Aortic Annuloplasty Device is indicated in adult patients undergoing bicuspid aortic valve repair for aortic insufficiency or concomitant to repair of an aortic aneurysm.

4. Annuloplasty Device Description

4.1. Overview

The HAART 200 Aortic Annuloplasty Device is a three dimensional annuloplasty ring designed to correct annular dilatation and/or maintain annular geometry of the aortic valve in adult patients undergoing bicuspid aortic valve repair. The Device consists of a titanium frame machined from medical grade Titanium 6Al-4V covered with medical grade polyester fabric affixed to the frame by suture.

4.2. Technological Characteristics

The HAART 200 Aortic Annuloplasty Device is comprised of three components: the implantable Annuloplasty Device, polyester Pledgets, and a Device Holder that is discarded during the procedure. Each of these components is briefly described below.

Annuloplasty Device

The HAART 200 Aortic Annuloplasty Device has a circular base geometry with two 10°, outwardly flaring, subcommissural posts positioned at 180° around the circumference of the device. The annuloplasty device consists of a titanium frame machined from medical grade Titanium 6Al-4V covered with medical grade polyester fabric affixed to the frame by suture. The Device materials and the manufacturing processes were specifically selected for use in an implantable medical device. The polyester fabric, ARF001, is manufactured for annuloplasty ring applications. HAART 200 Aortic Annuloplasty Devices are manufactured in 4 sizes ranging from 19mm to 25mm in 2mm increments.

The titanium frame of the Device provides the stiffness to reduce the diameter of the dilated aortic annulus while the polyester fabric provides material to support endothelialization and direct suturing of the subcommissural posts to the aortic valve annulus. The inner aspects of the Device posts have 2 layers of polyester fabric to facilitate suturing.

Polyester Pledgets

Polyester Pledgets are 7mm by 3mm and are made from the same ARF001 fabric used to cover the Device. The Pledgets are provided with the annuloplasty Device for use during surgery. They are supplied sterile in a separate package within the Device shelf box. Individual packages of 6 sterile pledgets are also available in the HAART Polyester Pledget 6-Pack.

Holder

The HAART 200 Aortic Annuloplasty Device is supplied on a Holder that may be attached to a Handle to facilitate positioning the Device during the procedure (Figure 2). The Holder is machined from polyphenylsulfone and is attached to the annuloplasty Device using a single suture (Figure 1). The Device can be removed from the Holder by severing the suture at any location along the face of the Holder.

Figure 1. HAART 200 Aortic Annuloplasty Device on holder


Figure 2. HAART 200 Aortic Annuloplasty Device and holder affixed to handle


Figure 3. HAART 201 Instrument Set

4.3. Accessories

The HAART 201 Instrument Set (Figure 3) is comprised of two (2) Handles, four (4) Sizers and a Gage Sphere. Handles are made from ASTM A276 stainless steel and the Sizers and Gage Sphere are fabricated from polyphenylsulfone.

The Handle may be threaded into the Sizers and Gage Sphere. The Handle may also be threaded into the face of the Holder to facilitate removal of the Device from the packaging and positioning the Device during the procedure (Figure 2). The Handle may be bent in the narrowed section to present the Sizers, Gage Sphere, and Device to the surgical site in the desired manner.

Warning: The Sizers, Gage Sphere, and Handle are intended for multiple uses provided they are inspected before each use for signs of damage. A cleaning and steam sterilization process has been validated for these reusable instruments for use by the hospital. Before each use, the Sizers and Gage Sphere should be visually inspected for crazing of polymer materials, cracks, signs of structural weakness, or unreadable

markings. The Handle should be inspected for visible cracks or signs of structural weakness before each use and after bending. Replace any instrument that exhibits these faults as they may not function properly and could cause patient injury. For more information on instruments, see the HAART 201 Instrument Set IFU.

5. Contraindications

- The Device is contraindicated in patients with a porcelain aorta.
- The Device is contraindicated in patients with evolving bacterial endocarditis.
- The Device is contraindicated in patients with heavily calcified valves.

6. Warnings

- The HAART 200 Aortic Annuloplasty Device is for Single Use Only. Do not re-use the Device. In addition to the risks listed under Complications, re-use may cause procedural complications including Device damage, compromised Device biocompatibility, and Device contamination. Re-use may result in infection, serious injury, or patient death.
- The decision to use an annuloplasty device must be made by the responsible physician on an individual basis after evaluation of the risks and benefits accrued to the patient in comparison to alternate treatment.
- Do not attempt to deform or reshape the annuloplasty Device.
- The HAART 200 Aortic Annuloplasty Device has been sterilized by gamma irradiation methods and is provided sterile in a double packaged container. No steam sterilization cycle has been validated for the sterilization of the Device.
- Valve repair patients subjected to subsequent dental procedures or other surgical procedures should receive prophylactic antibiotic drug therapy to minimize the risk of systemic bacteremia and prosthetic endocarditis.
- Correct annuloplasty device sizing is an important element of successful valve repair. Significant undersizing can result in valve stenosis or ring dehiscence. Oversizing can result in valve regurgitation. The size of the HAART 200 Aortic Annuloplasty Device is selected using Sizers consistent with the design intent of the Device. Use only the HAART Sizers included in the HAART 201 Instrument Set to select the proper Device size. Do not use the Holder as a sizing tool.

7. Precautions

- Only surgeons having received training in HAART Device implantation and sizing techniques should use this Device.
- Aortic valve leaflet reconstruction is routinely required to create valve competence during aortic valve repair. Annuloplasty using the HAART Device should be combined with leaflet reconstruction when indicated.
- Complex leaflet lesions including calcifications, fenestrations, perforations, nodular scarring and retraction, and tissue insufficiency requiring patch materials may be a risk factor for repair failure.
- Do not use the Device after the expiration date printed on the label.
- To ensure the sterility and integrity of the Device, it should be stored in the outer cardboard box until needed for introduction into the sterile field. Do not use a Device that has been removed from the double packaging and dropped, soiled, or otherwise damaged.
- Do not use the HAART 200 Aortic Annuloplasty Device if the tamper-evident seal is damaged, broken or missing.
- Dispose of used Devices as biohazardous waste.

- To avoid damage to the fabric covering the Device, do not use suture needles with cutting edges during implantation.
- Implantation sutures should be placed deeply to position the HAART devices at least 2 mm below the leaflet aortic junction and the tops of the subcommissural triangles in order to prevent contact between the leaflets and the Device.
- To prevent sutures from pulling through the native annular tissue, Pledgets should be used with the sutures as described in the Directions for Use. A sufficient number of broad horizontal mattress sutures should be used to eliminate gaps between the Device and the adjacent annular tissues and to support the tension associated with annular reduction.
- Sutures should be pulled tight and tied tightly so that the Device is in firm contact with the adjacent annular tissues. Suture tails should be tied down to the lateral aspect of the annular Pledgets, as described in the Directions for Use, to prevent leaflet damage from long suture tails. The sutures should also be cut very short.
- Provide careful monitoring of the patient's anticoagulation status when postoperative anticoagulation therapy is used. Surgeons who use the HAART 200 Aortic Annuloplasty Device should be current on anticoagulation regimens.

8. Magnetic Resonance (MR) Safety

MR Conditional

Non-clinical testing has demonstrated that the HAART 200 Aortic Annuloplasty Device is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3 T
- Maximum spatial field gradient of 4,000 G/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4.0 W/kg (First Level Controlled Operating Mode) at 3 T

RF Heating

Under the scan conditions defined above, the HAART 200 Aortic Annuloplasty Device is expected to produce a maximum temperature rise of less than 3.0°C after 15 minutes of continuous scanning.

Caution: The RF heating behavior does not scale with static field strength. Devices that do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artifact

In non-clinical testing, the image artifact caused by the device extends approximately 10 mm from the HAART 200 Aortic Annuloplasty Device when imaged with a gradient echo pulse sequence and a 3 T MRI system.

9. Potential Adverse Events

Each prospective patient should be informed about the benefits and risks of valve repair and annuloplasty surgery before the procedure. Serious complications, including death, are possible with any open heart surgery procedure including the implantation of the HAART 200 Aortic Annuloplasty Device. These potential complications include those associated with open heart surgery in general and the use of general anesthesia. The potential complications associated with the HAART 200 Aortic Annuloplasty Device and its implantation procedure are listed in Table 1.

Table 1. Possible Complications Associated with the Device/Procedure

Abrasion of the natural valve	Leaflet damage
Allergic reaction	Left Ventricular Outflow Tract Obstruction
Angina	Myocardial infarction
Aortic insufficiency	Neurological events (including TIA, stroke, and psychomotor deficit)
Arrhythmia	Pain (patient discomfort)
Atrioventricular block	Pericardial effusion
Death	Permanent pacemaker
Device explant	Pleural effusion
Device fracture	Psychological problems
Device migration or malposition requiring intervention	Renal insufficiency / failure
Endocarditis	Reoperation
Extended surgery time or aborted procedure	Respiratory system disorders
Fever	Ring dehiscence
Gastrointestinal disorders	Stenosis
Heart Failure	Suture injury to coronary arteries
Hematoma	Syncope
Hemolysis or hemolytic anemia	Thrombosis or thromboembolism
Hemorrhage	Toxic reaction
Hypertension	Wound healing problems
Infection – local, bacteremia, sepsis	

10. How Supplied

10.1. Packaging

The HAART 200 Aortic Annuloplasty Device is available in 19, 21, 23 and 25 mm sizes. Each HAART 200 Aortic Annuloplasty Device shelf box contains Pledgets and a single annuloplasty Device assembly consisting of the Device sutured onto its corresponding Holder. The assembled Device and Holder are packaged within nested, sealed trays. Pledgets are packaged separately within nested pouches. The packaging system is designed to ease placement of the Device into the sterile field. The components within the packaging are sterile if the pouches, trays, and lids are undamaged and unopened. The surfaces of the outer packaging are NONSTERILE and must not be placed in the sterile field.

10.2. Storage

Store the product in its original packaging, including the outer shelf box, in a clean and dry area to protect the product and minimize the potential for contamination. Stock rotation is recommended at regular intervals to ensure usage before the expiration date printed on the box label. Do not use the Device after the expiration date printed on the label.


11. Directions for Use

11.1. Sizing

The HAART 200 Aortic Annuloplasty Device size should be selected based on the free-edge length of the non-fused valve leaflet. The free-edge length is determined using the Sizers provided in the HAART 201 Instrument Set. Sizers are provided in 19, 21, 23, and 25 mm diameters to correspond to the four HAART 200 Aortic Annuloplasty Device sizes.

Correct annuloplasty device sizing is an important element of successful valve repair. The appropriate HAART 200 Aortic Annuloplasty Device size is selected by threading each individual Sizer onto the Handle and inserting it behind the non-fused valve leaflet such that the leaflet free-edge length between commissural insertions lies smoothly along the circumference of the Sizer. The appropriate Sizer for a given leaflet has been selected when the distance from one stippled area to the other matches the leaflet free-edge length from one commissure to the other (Figure 4). If the leaflet free-edge length is between two sizes, choose the smaller of the two size options. The size number on the corresponding Sizer indicates the appropriate Device size based on the leaflet measured.

Figure 4. Proper sizing of leaflet free-edge length using the sizer

To confirm the device size selected with the Sizer, the distance between the valve commissures can be measured using a valve sizer such as a Hegar dilator. The objective of this measurement is to ensure the diameter of the device chosen using the HAART 201 Sizers will not cause excessive reduction of the intercommissural distance.

11.2. Handling and Preparation Instructions

Each Device is supplied mounted on a Holder and is packaged in nested, sealed trays to ease transfer of the Device into the sterile field. Pledgets are packaged within nested pouches. The inner pouches and sealed trays should be inspected for damage prior to opening. Do not use the Device if the sterile packaging has been compromised.

The Holder should be attached to the Handle contained within the HAART 201 Instrument Set to facilitate removal of the Device from the packaging and placement of the Device within the aortic root.

11.3. Device Implantation


Post Sutures

Insertion is begun by suturing the two posts of the Device to the two subcommissural areas using “Cabrol-like” configurations with generous

bites taken in the aortic wall, using 4-0 or 3-0 coated braided polyester horizontal mattress sutures supported by Pledgets above the annulus (Figure 5).


The post sutures are placed with the HAART 200 Aortic Annuloplasty Device on the Holder above the valve (Figure 6). The post sutures should be placed so that the tip of the device post will be positioned at least 2mm below the top of the subcommissural space.

Figure 5. Suturing technique for the HAART 200 Aortic Annuloplasty Device Posts

The HAART 200 Aortic Annuloplasty Device is covered with polyester fabric, which allows endothelialization. Horizontal mattress post sutures are used primarily to correctly position the Device prior to placement of additional looping sutures around the leaflet sections

Figure 6. Placing the post sutures

of the Device. Only the horizontal mattress post sutures should catch the fabric and then only just on the inside of the posts (Figure 5). Extra fabric has been added in those areas to allow needle passage, but only very superficial bites are required.

After the two post sutures are placed, the Device is lowered below the native valve, and the suture attaching the Device to the Holder is cut. The HAART 200 Aortic Annuloplasty Device must be removed gently from the Holder, primarily by pushing the Device off the Holder near the center of each leaflet section of the Device. Only after the Device has been pushed off the Holder should the Holder be removed from beneath the valve.

Sinus Section Sutures

Sinus section sutures loop around the Device and pass up through the annulus, taking deep bites into the aorta and emerging above the valve onto fine Pledgets (Figures 7 and 8). Three looping sutures are placed around the non-fused sinus segment of the Device and four (4) looping sutures are placed around the sinus segment corresponding to the fused leaflet, two (2) on either side of the raphe (Figure 9).

Suture Management

After all nine sutures are placed, each is tied firmly over the Pledget with 8 knots, ensuring that the leaflet posts are buried back into the subcommissural regions and the Device fabric is kept below leaflet tissues. The knots should be thoroughly tightened to prevent the sutures from coming untied.

Figure 7. Sinus section looping sutures


Figure 8. Suturing of the sinus sections with the device below the valve leaflets.



Figure 9. Final suture configuration


Long annular suture tails in the coronary sinus can cause leaflet injuries. As a final step of annular suturing, one or both needles from the tied annular suture should be passed downward through the center of the lateral aspect of the Pledget, and the suture should be again tied down to the Pledget with 6 more knots. This maneuver directs the suture tails down and away from the leaflets.

Figure 10 illustrates the procedure for proper management of the suture tails.

Panel A: The horizontal mattress sutures that hold the Device firmly up under the annulus emerge above the annulus and are supported by fine polyester Pledgets. The sutures are tied tightly with 8 knots over the Pledgets (arrow). These suture lines are reducing annular size significantly, and therefore, can be associated with significant tension. Thus, extra care should be taken with tying good 8-throw knots that will not come untied.

Figure 10. Management of Annular Sutures

Panel B: Each suture is passed down through the lateral Pledget (arrow) and tied again with 6 knots. This step positions the final knot and suture tails laterally and under the Pledget.

Panel C: The suture tails are cut very short beneath the Pledgets (arrow).

Panel D: The double knot prevents each suture from coming untied, and with the second knot beneath the Pledget, the suture tails are directed laterally and down into the annulus, preventing contact with the leaflets.

At the end of the Device implant procedure, all annular sutures should be carefully inspected by the surgeon, and if any question exists about the position of a given suture tail, the knot tower should again be sutured down and away from the leaflet with a fine 6-0 polypropylene suture.

Leaflet Reconstruction

When indicated, leaflet reconstruction should be performed after annuloplasty Device insertion and according to established aortic leaflet reconstruction techniques. Annuloplasty device insertion moves the annulus and leaflets toward the center of the valve, but significant leaflet prolapse may still be present. Leaflet prolapse may be corrected by plication sutures in the leaflet free-edge, or by closing the cleft of a Type 1 valve. If present,

structural leaflet defects may be corrected, using pericardial reconstruction or other methods, consistent with the surgeon's training and preferences. The goal is to achieve equal lengths of both leaflets and equal effective heights of > 8 mm.

Settepani F, Cappai A, Raffa G, et al. Cusp repair during aortic valve-sparing operation: technical aspects and impact on results. *J Cardiovasc Med (Hagerstown)*. 2015 Apr;16(4):310-7.

Mazzitelli D, Stamm C, Rankin JS, et al. Leaflet reconstructive techniques for aortic valve repair. *Ann Thorac Surg*. 2014 Dec;98(6):2053-60.

Assessment of Valve Repair

After completing the Device implantation and leaflet reconstruction, the leaflets should be inspected and the aorta at the commissures should be stretched laterally to assess relative leaflet lengths and heights. Both leaflets should be equal in length and vertical with good effective height and surface of coaptation. The leaflets should meet in the midline with no evidence of central gaps or leaflet prolapse. Finally, leaflets should be opened widely to ensure a good orifice.

Effective height and geometric height may be assessed using the Gage Sphere provided in the HAART 201 Instrument Set (Figure 11). The Gage Sphere has two different vertical scales marked on the surface. One scale originates at the apex of the Sphere and is intended for estimating the vertical effective height of the valve leaflet from the leaflet base to the free-edge margin (Figure 12). For a successful repair, leaflet effective height should be approximately 8 to 10 mm. Thus, as the Gage Sphere is gently pressed down into the leaflet-sinus complex, a successful repair will be associated with the leaflet free-edge being at the level of the Gage Sphere equator. The second scale may be used to assess leaflet geometric lengths, when desired, for further evaluation of leaflet size and symmetry after the overall valve repair.

11.4. Sterilization

The HAART 200 Aortic Annuloplasty Device is provided sterile on the Holder and must not be resterilized. Devices that have been damaged or contaminated should not be used. Pledgets are provided sterile and must not be resterilized. Pledgets that have been damaged or contaminated from patient contact should not be used.

Figure 11. Diagram of leaflet geometric height and effective height



Figure 12. Geometric Height Scale (left) and Effective Height Scale (right)

Disclaimer of Warranties

Although the HAART 200 Aortic Annuloplasty Device and HAART 201 Instrument Set, hereafter referred to as "product," have been manufactured under carefully controlled conditions, Corcym has no control over the conditions under which this product is used. Corcym and its affiliates (collectively, "Corcym"), therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Corcym shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Corcym to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty and Limitation of Liability is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty and Limitation of Liability shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty and Limitation of Liability did not contain the particular part or term held to be invalid.

Patents

Patents: US8,163,011; US8,425,594; US9,161,835; US9,814,574; US9,844,434; US10,130,462; US10,327,891; CA 2,665,626; JP5881653; JP5877205; JP6006218; EP2621407; EP2621408; other applications pending.

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Corcym Inc.
Corcym S.r.l.