Instructions for Use - HAART 201 Instrument Set

Instructions for Use -

HAART 201 Instrument Set2

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Instructions for Use - HAART 201 Instrument Set

Table of Contents

1. Instrument Labeling Symbols	2
2. Indications	2
3. Description	2
4. Warnings	4
5. Precautions	4
6. How Supplied	4
7. Directions for Use	5
8. Sizer and Gage Sphere Accuracy Statement	5
9. Cleaning and Sterilization Prior to Use	5
Disclaimer of Warranties	7
Patents	8
Figures.....	137

1. Instrument Labeling Symbols

Manufacturer

Nonsterile

Batch code

Catalog number

Authorized European Rep

Product complies with requirements of directive 93/42/EEC for medical devices

Consult Instructions for Use

2. Indications

The HAART 201 Instrument Set is intended to be used as an aid in the implantation of the HAART 200 Aortic Annuloplasty Device.

3. Description

The HAART 201 Instrument Set consists of four (4) Sizers, one (1) Gage Sphere, and two (2) Handles provided in a stainless steel tray (Figure 1 on page “Obrázok 1 Súprava nástrojov HAART 201” on page 137). The Instrument Set is used to select the proper HAART 200 Aortic Annuloplasty Device size for a given patient and to assist

in the assessment of leaflet coaptation after completion of the overall repair of the aortic valve. The HAART 201 Instrument Set is reusable; however, it must be cleaned and sterilized by autoclave (steam) prior to each use. Refer to the sections below for detailed information on accessory use and sterilization.

3.1 Handle

The Handle is machined from ASTM A276 stainless steel, and cleaned and passivated per ASTM specification F86. The Handle may be threaded into the Sizers, Gage Sphere, and Holder to facilitate positioning during the procedure (Figure 2 on page “Obrázek 2 Zařízení k aortální anuloplastice HAART 200 a držák připojené k rukojeti” on page 137). The stainless steel Handle may be bent in the narrowed section to present the Sizers, Gage Sphere, and Device to the surgical site in the desired manner. The Handle is very similar in design and manufacture to handles supplied with replacement heart valve prostheses and is intended to be used in a similar manner.

3.2 Sizers

The Sizers are machined from polyphenylsulfone. Each Sphere has a through hole along its axis (Figure 3 on page “Obrázek 3 Ilustrace měřicí pomůcky” on page 138). The top of the through hole is threaded to receive the instrument Handle.

The Sizers are laser marked with a line at the equator and 15° wide stippled regions, oriented 180° apart, to assist the user with assessment of leaflet free-edge length. The size of the corresponding HAART 200 device is prominently marked on the upper surface of each Sizer.

3.3 Gage Sphere

The Gage Sphere is machined from polyphenylsulfone and has a through hole along its axis. The top of the through hole is threaded to receive the instrument Handle.

After completion of leaflet repair, effective height and leaflet geometric height may be assessed using the Gage Sphere provided in the HAART 201 Instrument Set. The Gage Sphere has two different vertical scales marked on the surface (Figure 4 on page “Obrázek 4 Škála k měření délky cípů (vlevo) a škála k měření efektivní výšky (vpravo)” on page 139). One scale originates at the apex of the Sphere and is intended for estimating the effective height of the valve leaflet. For a successful repair, leaflet effective height should be approximately 8 to 10 mm. The second scale may be used to assess leaflet geometric heights, when desired,

for further evaluation of leaflet size and symmetry after the overall valve repair.

4. Warnings

- Do not use other manufacturer's annuloplasty sizers. Other annuloplasty sizers may not indicate the appropriate HAART 200 Aortic Annuloplasty Device size.
- Correct annuloplasty device sizing is an important element of successful valve repair. Significant undersizing can result in valve stenosis or ring dehiscence. Oversizing can result in valve regurgitation. The size of the HAART 200 Aortic Annuloplasty Device is selected using Sizers consistent with the design intent of the Device. Use only the HAART Sizers included in the HAART 201 Instrument Set to select the proper Device size.
- Use only the HAART Handle to interface with the Device Holder, Sizers and Gage Sphere.
- The Sizers, Gage Sphere, and Handle are provided NONSTERILE and must be cleaned and sterilized prior to each use.
- Only surgeons having received training in HAART device implantation and sizing techniques should use this Device.
- The HAART 201 Instrument Set is not intended for use with devices other than the HAART 200 Aortic Annuloplasty Device.
- Ultrasonic cleaning should not be used on the HAART 201 Sizers or Gage Sphere since it may cause crazing of the polymer material and markings.

5. Precautions

The Sizers and Gage Sphere should be visually inspected for crazing of polymer materials, cracks and signs of structural weakness, or unreadable markings before each use.

The Handle should be inspected for visible cracks or signs of structural weakness before each use and after bending. Replace any instrument that exhibits these faults as they may not function properly and could cause patient injury.

6. How Supplied

The HAART 201 Instrument Set is provided NONSTERILE in a stainless steel tray. The instruments are provided clean, but NONSTERILE.

Upon initial receipt of the Instrument Set, remove all packaging materials and carefully examine all instruments for signs of shipping damage. The instruments must be visually inspected, cleaned, and sterilized before each use.

Discard any instruments that show signs of damage or unreadable markings.

Store the instrument set in a clean environment.

7. Directions for Use

The HAART 201 Instrument Set is intended to be used as an aid in the implantation of the HAART 200 Aortic Annuloplasty Device. See the Directions for Use for the HAART 200 Aortic Annuloplasty Device.

8. Sizer and Gage Sphere Accuracy Statement

Sizers serve as a visual gage for choosing an appropriate Device size based on the leaflet free-edge lengths of the aortic valve. Size diameters are accurate to ± 0.010 inches (0.254 mm) of the specified size.

The Gage Sphere is intended only as a reference and visual aid to assess leaflet effective height and leaflet geometric length following leaflet repair. The scales on the Gage Sphere are accurate to ± 0.008 inches (0.2 mm) between adjacent marks.

9. Cleaning and Sterilization Prior to Use

Warning: The HAART 201 Instrument Set is provided NONSTERILE and must be thoroughly cleaned and sterilized prior to use.

9.1 Cleaning Instructions

The Handle should be removed from the Sizer or the Gage Sphere prior to cleaning.

Note: If transport of the instruments to the processing area is likely to be delayed, place the instruments in a covered container with an appropriate detergent (e.g. Enzol Enzymatic Detergent or equivalent) to delay drying.

Within a maximum of 2 hours after use, remove gross soil using a cold (25° to 30°C) water rinse with running tap water and soak for 1 minute in an enzyme detergent (e.g. Enzol Enzymatic Detergent) prepared with 2 oz. per gallon of water. Immerse instruments in a container of an appropriate cleaner (e.g. Enzol) prepared with 2 ounces to 1 gallon tap water (15° to 30°C) for a minimum of two (2) minutes. Using a bristle brush, clean the threaded area of the sizers and gage spheres for 15 to 20 seconds to ensure removal of any impacted soil.

Rinse in running tap water (20° to 35°C) until all traces of cleaning solution are visually removed. Visually inspect for any remaining soil and repeat the process if necessary. Dry using a soft cloth and transport to the cleaning step.

Clean the instrument in an automated washer using the following parameters:

Automated Cleaning Parameters

Process Step	Time (mm:ss)	Temp (C)	Cleaner
Enzymatic Wash	04:00	Hot tap water	Steris Prolystica 2X Concentrate Enzymatic Presoak and Cleaner 3 mL/liter
Wash	02:00	60	Steris Prolystica Ultra-Concentrate Neutral Detergent 0.7 mL/liter
Rinse	02:00	70	N/A
Dry	15:00	Low	N/A

Visually inspect each instrument for remaining soil and moisture. If soil remains, repeat the cleaning process. If instruments are wet, use filtered compressed air or a lint free wipe to dry.

Warning: Ultrasonic cleaning should not be used on the Sizers since it may cause crazing of the polymer material and markings.

The user is responsible for the qualification of any deviation from the recommended cleaning method.

9.2 Sterilization Instructions

The Handles should be removed from the Sizer or Gage Sphere prior to cleaning and sterilization.

The four Sizers, the single Gage Sphere, and Handles should be placed in the instrument tray for steam sterilization. The following sterilization cycle is recommended:

Temperature: 132°C (270°F)

Dwell time: 4 minutes

Cycle: Pre-vacuum

Dry time: 20 minutes minimum

The user is responsible for qualification of any deviations from the recommended sterilization method.

Corcym acknowledges that the sterilization cycle described below is routinely used by European hospitals and has been validated for many similar heart valve sizer systems and surgical instruments.

Pre-vacuum: cycle

Temperature: 134°C (274°F)

Dwell time: 3 minutes

Dry time: 20 minutes

The components of the HAART 201 Instrument Set do not pose any significant challenge to sterilization using this widely validated process. The instrument materials can withstand temperatures up to 137°C.

9.3 Re-use

Corcym has performed a minimum of 50 sterilization and use cycles on the HAART 201 Instrument Set without any signs of wear being noted.

Caution: The Sizers, Gage Sphere, and Handle are intended for multiple uses as long as they are inspected before each use for signs of damage. The Sizers and Gage Sphere should be visually inspected for crazing of polymer materials, cracks and signs of structural weakness or unreadable markings before each use. The Handle should be inspected for visible cracks or signs of structural weakness before each use and after bending. Replace any instrument that exhibits these faults, as it may not function properly and could cause patient injury.

Disclaimer of Warranties

Although the HAART 200 Aortic Annuloplasty Device and HAART 201 Instrument Set, hereafter referred to as “product,” have been manufactured under carefully controlled conditions, Corcym has no control over the conditions under which this product is used. Corcym and its affiliates (collectively, “Corcym”), therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Corcym shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Corcym to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty and Limitation of Liability is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty and Limitation of Liability shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty and Limitation of Liability did not contain the particular part or term held to be invalid.

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Corcym Inc.
Corcym S.r.l.